pH - Keep it Wet!!! - a pH FAQ

Frequently Asked Questions

Q: How often do I need to calibrate my pH tester or meter?

A: Depending on the accuracy required, once a day if several measurements are required during the day; or at least once weekly if you use 2 or 3 times per week.

Q: Why does the pH tester or meter Display "1" or show erratic reading when power on?

A: When the electrode is held in the air and not dipped or in contact with a solution, it will display erratic readings. This is because the reference and glass electrode are in close contact, therefore the electronics cannot show a fixed reading. Dip the electrode into any solution and a reading will appear.

Q: Why does the pH tester or meter take sometime to get a stable reading?

A: Electrode was dry or not properly conditioned. Reference junction is possibly clogged up. Remedy: Pre-soak the electrode for 1 hour in tap water before using.

Q: How long will a pH tester last?

A: All pH glass electrode will age with time and temperature.

Most glass electrodes have a storage life of about 1 year. A life-span of about 3 to 6 months for frequent/continuous use.

Q: What other problems are associated with measurement?

A: Temperature is important as liquid properties are affected. Therefore, a correct temperature of 25 °C for Buffer Solutions should be maintained during calibration. However correction of values may be made if the temperature of the solution is known, by comparing with the table below. Different temperatures of test samples will give inconsistent results.

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