



# Non-contact FOREHEAD THERMOMETER

**No-touch infrared technology**  
For safety and hygiene.

**Forehead use**  
Displays "oral equivalent"  
temperature.

**Simple to use**  
Point the thermometer at  
the forehead and press START  
It's that easy!

**High temperature alert**  
Audible alert when  
temperatures  
are  $> 37.5^{\circ}\text{C}$  ( $99.5^{\circ}\text{F}$ )

**Memory**  
Stores the 25 most recent  
readings for easy recall.

**FDA-cleared**  
For use in adults and children.

**Quick one second readings**  
Instant results - no waiting.  
Beeps when done.

**Mute setting**  
Disable audible alert  
that could disturb a  
sleeping child.

**Bright LED backlight**  
Makes temperatures  
easy to read, even in  
the dark.



## NON-CONTACT, ONE-SECOND READINGS FOR QUICK, SAFE SCREENING.

\*FDA-cleared for safety & accuracy for kids and adults. ETI non-contact forehead thermometer measures infrared energy radiated from the skin at the centre of the forehead area. This captured energy is collected through the lens and converted to a body temperature value (displays oral equivalent temperature).

It's helpful to know each individual's normal temperature when they are well. This is the best way to accurately diagnose a 'high temperature'. To get a reference value, we suggest taking multiple readings (5 to 10) over one minute, taking the highest number to determine "normal temperature".